CMV-IgM-ELA Test PKS medac

Deutsch/English/Français

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CMV-IgM-ELA Test PKS medac

Enzyme immunoassay with Pipetting Control System (PCS) for the detection of IgM antibodies to Cytomegalovirus (CMV)

Cat. no.: 110-PKS

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION


Cytomegalovirus (CMV) belongs to the family of human pathogenic Herpesviridae which is characterised by a double-stranded DNA genome. It is characteristic for these viruses to persist latently in the organism after primary infection. Reactivations of the virus can therefore occur under certain circumstances.

The course, CMV infections take in immunocompetent individuals, is normally without any symptoms. In individuals who suffer from immunosuppression (e.g. transplantation patients, HIV infected individuals, tumor patients, newborns) a variety of serious symptoms are observed.


CMV is one of the most frequent pathogens in prenatal infection which can cause a variety of serious disorders, also including late damages in children born without symptoms.

The CMV-IgM-ELA test PCS medac can be applied for the detection of specific CMV IgM antibodies in serum samples. An acute CMV infection can be diagnosed by the demonstration of IgM antibodies. It has to be taken into consideration that IgM antibodies are not always produced. In immunosuppressive patients, for example, CMV IgM antibodies can rarely be detected during reactivations and reinfections.


TEST PRINCIPLE


The plate is coated with anti-human IgM immunoglobulin.

IgM of the specimen is selectively bound to the wells.

The CMV-specific IgM antibodies bind to the peroxidase-labelled CMV antigen (AG = antigen, P = peroxidase).

Incubation with TMB-substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ No unspecific reactions and no false positive results caused by rheumatoid factors.
- ☞ No blocking of IgM antibodies by high IgG titer.
- ☞ The Pipetting Control System allows to monitor visually each pipetting step through colour changes.
- ☞ The breakable microwell strips permit efficient use of the test.
- ☞ Suitable for automation on open ELISA devices.

KIT CONTENTS

Cat. no.: 110-PKS

1.

| |
|------------|
| MTP |
|------------|

Microplate: 12 x 8 wells (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with goat anti-human IgM immunoglobulin, BSA and pH indicator, ready to use.
2.

| | |
|----------------|----------|
| CONTROL | - |
|----------------|----------|

Negative control: 1 vial with 1.5 ml, human serum, ready to use, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
3.

| | |
|----------------|----------|
| CONTROL | + |
|----------------|----------|

Positive control: 1 vial with 1.5 ml, human serum, ready to use, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
4.

| |
|-----------|
| WB |
|-----------|

Wash buffer: 1 bottle with 100 ml, PBS/Tween (10x), pH 7.2 - 7.4, contains ProClin™ 300.
5.

| |
|----------------|
| VIR-DIL |
|----------------|

Sample diluent: 1 bottle with 110 ml, PBS/Tween/BSA, pH 7.2 - 7.4, ready to use, contains ProClin™ 300.
6.

| |
|----------------|
| ANTIGEN |
|----------------|

CMV-IgM-ELA (Enzyme Labelled Antigen): 2 vials with 5.0 ml each, lyophilised, contains FCS, stained red, HRP-conjugated.
7.

| |
|------------|
| TMB |
|------------|

TMB-substrate: 1 vial with 10 ml, ready to use.
8.

| |
|-------------|
| STOP |
|-------------|

Stop solution: 2 vials with 11 ml each, 0.5 M sulfuric acid (H₂SO₄), ready to use.

1. STORAGE AND STABILITY

| MATERIAL/REAGENT | STATE | STORAGE | STABILITY |
|------------------|---------------|-------------------------------|-------------------|
| Test kit | unopened | 2...8°C | until expiry date |
| Microplate | opened | 2...8°C in bag with desiccant | 6 weeks |
| Controls | opened | 2...8°C | 6 weeks |
| Wash buffer | diluted | 2...8°C | 6 weeks |
| Sample diluent | opened | 2...8°C | 6 weeks |
| CMV-IgM-ELA | reconstituted | 2...8°C | 5 days |
| | | ≤ -18°C * | 6 weeks |
| TMB-substrate | opened | 2...8°C | 6 weeks |
| Stop solution | opened | 2...8°C | until expiry date |

* Aliquote and do not freeze again after thawed once!

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature to prevent condensation.

Calculate the number of wells required.

3.1. Microplate

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated under point 1.

Note: The microplate wells have a light green colour. Eventually occurring greenish brown stains inside the wells are due to the production process and do not influence the test performance.

3.2. Wash buffer

Mix one volume of wash buffer (10x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

3.3. CMV-IgM-ELA

Reconstitute the lyophilised CMV-IgM-ELA with 5.0 ml sample diluent each. Mix gently and take care that particles sticking to the closure are also dissolved.

After reconstitution, the CMV-IgM-ELA has a red colour and is ready to use.

Do not mix test specific reagents (microplate, controls, CMV-IgM-ELA) from different kit lots.

In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all virological ELISA medac.

Reagents from other manufacturers must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed .

4. SPECIMEN

4.1. The test is suitable for serum samples.

4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.

4.3. Sera have to be diluted 1 : 100 with sample diluent. They can be diluted further for titer determination.

5.A. TEST PROCEDURE

5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready to use and do not have to be pre-washed.

- 5.2. Leave well A1 empty as blank (see 6.A.). Add 50 µl each of the diluted sample as well as negative control and positive control in duplicates to the wells.

After pipetting the samples (pH neutral or basic fluid) the wells turn blue. A missing colour change in one well indicates that no sample or control has been added.

If necessary the microplate wells can be kept at RT in a humid chamber up to 30 min or at 2 - 8 °C up to 60 min before proceeding.

- 5.3. Incubate the microplate wells for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. Reconstitute the CMV-IgM-ELA (see 3.3.).
- 5.5. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

- 5.6. Add CMV-IgM-ELA (coloured red) to each well (except A1).

50 µl of ELA have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated devices, 60 µl of ELA have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.

The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.

- 5.7. Incubate again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.8. After incubation wash microplate wells again (see 5.5.).
- 5.9. Add 50 µl of TMB-substrate to each well (also A1) and incubate for 30 min (\pm 2 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.
- 5.10. Stop the reaction by adding 100 µl of stop solution to each well (also A1). Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the stop solution!

5.B. TABLE FOR THE TEST PROCEDURE

| | Blank (A1) | Negative control | Positive control | Sample |
|---|------------|------------------|------------------|-------------|
| Negative control | - | 50 µl | - | - |
| Positive control | - | - | 50 µl | - |
| Sample | - | - | - | 50 µl |
| Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer. | | | | |
| CMV-IgM-ELA | - | 50/60 µl *) | 50/60 µl *) | 50/60 µl *) |
| Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer. | | | | |
| TMB-substrate | 50 µl | 50 µl | 50 µl | 50 µl |
| Incubate for 30 min at 37 °C in the dark. | | | | |
| Stop solution | 100 µl | 100 µl | 100 µl | 100 µl |
| Photometric reading at 450 nm (ref. 620 - 650 nm) | | | | |

*) manual/automatic procedure (see 5.6.)

6.A. CALCULATION OF RESULTS (VALIDITY)

- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.
- * The mean OD value of the **negative control** has to be < 0.100.
The mean OD value of the **positive control** has to be > 0.800.
- * **Cut-off = mean OD value of the negative control + 0.140**
- * **Grey zone = Cut-off ± 10 %**

Repeat the run if the results do not meet the specification!

6.B. INTERPRETATION OF RESULTS

* Samples with OD values below the lower limit of the grey zone are reported as **NEGATIVE**.

* Samples with OD values within the grey zone are reported as **EQUIVOCAL**.

These should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.

* Samples with OD values exceeding the upper limit of the grey zone are reported as **POSITIVE**.

* The results should always be interpreted in connection with clinical data and additional diagnostic parameters.

* Cross reactions, caused by antibodies against other Herpesviridae, cannot be excluded in single cases.

* Very high concentrations of lipids in negative samples can cause false positive results.

High concentrations of Hemoglobin or Bilirubin do not influence the test results.

7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation.

7.A. SENSITIVITY AND SPECIFICITY

482 serum specimen, selected from blood donors and patients, were tested with the CMV-IgM-ELA test PCS medac in comparison to a highly specific ELISA that is used routinely in laboratories. The results obtained are shown in the table below.

| CMV-IgM-ELA Test PKS | Reference test | |
|----------------------|----------------|----------|
| | negative | positive |
| | negative | 314 |
| positive | 1 | 166 |

Sensitivity = 99.40 %

Specificity = 99.68 %

Positive predictive value: 99.40 %

Negative predictive value: 99.68 %

7.B. PRECISION

| Sample | Intra-assay variation | | | | Sample | Inter-assay variation (n = 11) | | |
|--------|-----------------------|-------|--------|----|--------|-----------------------------------|-------|--------|
| | mean OD | SD | CV (%) | n | | mean OD | SD | CV (%) |
| NC | 0.034 | 0.005 | 15 | 24 | NC | 0.050 | 0.010 | 19 |
| BC | 0.314 | 0.013 | 4 | 24 | BC | 0.230 | 0.021 | 9 |
| PC | 2.269 | 0.048 | 2 | 24 | PC | 1.456 | 0.064 | 4 |
| N° 1 | 0.049 | 0.008 | 15 | 27 | N° 4 | 0.058 | 0.007 | 13 |
| N° 2 | 0.317 | 0.023 | 7 | 27 | N° 5 | 0.264 | 0.037 | 14 |
| N° 3 | 2.673 | 0.086 | 3 | 27 | N° 6 | 0.271 | 0.039 | 14 |
| | | | | | N° 7 | 1.788 | 0.073 | 4 |

NC = negative control; BC = weak positive control (not included in the kit);
PC = positive control

GENERAL HANDLING ADVICES

- * To avoid cross contamination do not exchange the vials and their screw caps.
- * The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.
- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents) should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Flik, J.: A comparison of 2 quantitative ELISAs used to detect IgG antibodies against HCMV in transplant recipients. Poster presented at the 5th Intern. Cytomegalovirus Conference, (1995).

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Steinmann, J. and Bischoff, J.: Comparison of serological methods for the detection of *Cytomegalovirus* infection. Lab. Med. 15, 585 - 589, (1991).

Tönnies, R., Flik, J., Franke, D., Metzger, C., Daiminger, A., Bäder, U. and Enders, G.: Comparison of three methods used to differentiate between primary and recurrent or long-term human *Cytomegalovirus* infections in pregnant women. Poster presented at the 6th International Cytomegalovirus Workshop, Orange Beach, AL, USA, (1997).

Weber, B., Prosser, F., Munkwitz, A. and Doerr, H. W.: Serological diagnosis of *Cytomegalovirus* infection: comparison of 8 enzyme immunoassays for the detection of HCMV-specific IgM antibody. Clin. Diagn. Virol. 2, 245 - 259, (1994).