## Rubella-IgM-ELA Test PKS medac

---

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## **LITERATURHINWEISE**

s. S. 22

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<b>Rubella-IgM-ELA Test PKS medac</b>
Enzyme immunoassay with Pipetting Control System (PKS) for the detection of IgM antibodies to rubella virus

Cat. no.: 135-PKS

FOR IN VITRO DIAGNOSTIC USE ONLY

### **INTRODUCTION**

The pathogen causing rubella (German measles), rubella virus, belongs to the family of *Togaviridae* and is assigned to the genus *Rubivirus*. The rubella virus has a single-stranded RNA genome and occurs only in one serotype.


Rubella appears mainly in childhood and takes normally a mild course. The disease is characterised by a small-spotted exanthema and is accompanied by lymphadenosis.

During the first trimenon of pregnancy, infections with the rubella virus can cause a foetal infection, which can lead to irreversible embryofoetal damages.


Rubella virus-specific IgM antibodies can be detected 4 to 30 days after eruption of the rash and persist for 1 to 4 months. Therefore the IgM determination with the Rubella-IgM-ELA test PCS medac is a suitable method for diagnosis of an acute rubella virus infection.

Occasionally, IgM antibodies can persist up to one year or do not appear at all. Therefore the IgM result should always be interpreted together with anamnestic data and possible further examinations of the rubella virus infection.


## TEST PRINCIPLE


The plate is coated with anti-human IgM immunoglobulin.

IgM of the specimen is selectively bound to the wells.

The complex of mAb-peroxidase-conjugate and rubella virus antigen binds to rubella virus-specific IgM antibodies (AG = antigen, P = peroxidase).

Incubation with TMB-substrate (\*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

### Advantages of the test

- ☞ No unspecific reactions and no false positive results caused by rheumatoid factors.
- ☞ No blocking of IgM antibodies by high IgG titer.
- ☞ The Pipetting Control System allows to monitor visually each pipetting step through colour changes.
- ☞ The breakable microwell strips permit efficient use of the test.
- ☞ Suitable for automation on open ELISA devices.

## **KIT CONTENTS**

**Cat. no.: 135-PKS**

1. **MTP**  
Microplate: 12 x 8 wells (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with goat anti-human IgM immunoglobulin, BSA and pH indicator, ready to use.
2. **CONTROL -**  
Negative control: 2 vials with 0.75 ml each, human serum, ready to use, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
3. **CONTROL +**  
Positive control: 2 vials with 0.75 ml each, human serum, ready to use, contains, BSA, phenol, ProClin™ 300 and gentamycine sulfate.
4. **WB**  
Wash buffer: 1 bottle with 100 ml, PBS/Tween (10x), pH 7.2 - 7.4 contains ProClin™ 300.
5. **VIR-DIL**  
Sample diluent: 1 bottle with 110 ml, PBS/Tween/BSA, pH 7.2 - 7.4, ready to use, contains ProClin™ 300.
6. **ANTIGEN-DIL**  
Antigen diluent: 1 vial with 14 ml Tris/NaCl, ready to use, contains ProClin™ 300.
7. **ANTIGEN**  
Antigen: 4 vials with 3.0 ml each, rubella virus antigen/BSA, lyophilised.
8. **CON**  
Conjugate: 1 vial with 0.3 ml, anti-rubella virus IgG antibody, monoclonal (mouse), HRP-conjugated, stained red, ready to use, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
9. **TMB**  
TMB-substrate: 1 vial with 10 ml, ready to use.
10. **STOP**  
Stop solution: 2 vials with 14 ml each, 0.5 M sulfuric acid (H<sub>2</sub>SO<sub>4</sub>), ready to use.

## **1. STORAGE AND STABILITY**

<b>MATERIAL/REAGENT</b>	<b>STATE</b>	<b>STORAGE</b>	<b>STABILITY</b>
Test kit	unopened	2...8°C	until expiry date
Microplate	opened	2...8°C in bag with desiccant	6 weeks
Controls	opened	2...8°C	6 weeks
Wash buffer	diluted	2...8°C	6 weeks
Sample diluent	opened	2...8°C	6 weeks
Antigen diluent	opened	2...8°C	6 weeks
Antigen-conjugate mixture	ready to use	2...8°C	1 week
TMB-substrate	opened	2...8°C	6 weeks
Stop solution	opened	2...8°C	until expiry date

Do not use the reagents after the expiry date.

## **2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED**

- 2.1. Water for injection (H<sub>2</sub>O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of Wash Buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

## **3. PREPARATION OF THE REAGENTS**

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

### **3.1. Microplate**

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated in table 1.

**Note: The microplate wells have a light green colour. Eventually occurring greenish brown stains inside the wells are due to the production process and do not influence the test performance.**

### 3.2. Wash buffer

Mix one volume of wash buffer (10 x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10 x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

**Crystals in the wash buffer (10 x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.**

### 3.3. Antigen-conjugate mixture

Reconstitute the lyophilised antigen with 3.0 ml **antigen diluent** each. Mix gently and take care that particles sticking to the closure are also dissolved. Add 50 µl conjugate to 3.0 ml antigen 60 min before use.

**After reconstitution of the antigen and addition of the conjugate, the antigen-conjugate mixture has a red colour and is ready to use.**

**The ready-to-use antigen-conjugate mixture can be stored at 2 - 8 °C for one week (see 1.)! Attention! Do not freeze!**

**Do not mix test specific reagents (microplate, controls, conjugate, antigen, antigen diluent) from different lots.**

**In contrast to that sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all virological ELISA medac.**

**Reagents from other manufacturers must not be used in general.**

**Valid and reproducible results are only obtained if the test procedure is precisely followed .**

## 4. SPECIMEN

4.1. The test is suitable for serum samples.

4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.

4.3. Sera have to be diluted 1 : 100 with sample diluent. They can be diluted further for titer determination.

## **5.A. TEST PROCEDURE**

- 5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

**Microplate wells are ready to use and do not have to be pre-washed.**

- 5.2. Leave well A1 empty as blank (see 6.A.). Add 50 µl each of the diluted sample as well as negative control and positive control in duplicates to the wells.

**After pipetting the samples (pH neutral or basic fluid) the wells turn blue. A missing colour change in one well indicates that no sample or control has been added.**

**If necessary, the microplate wells can be kept in a humid chamber up to 60 min at 2 - 8 °C before proceeding.**

- 5.3. Incubate the microplate wells for 60 min ( $\pm$  5 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.4. Prepare the antigen-conjugate mixture (see 3.3.).
- 5.5. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

**Do not allow the wells to dry out! Proceed immediately!**

- 5.6. Add antigen-conjugate mixture (coloured red) to each well (except A1).

**50 µl of antigen-conjugate mixture have to be pipetted into the wells if the test is done manually.**

### **Please note:**

**When working with automated devices, 60 µl of antigen-conjugate mixture have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.**

**The suitability of the test automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.**

- 5.7. Incubate again for 60 min ( $\pm$  5 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil.
- 5.8. After incubation wash microplate wells again (see 5.5.).
- 5.9. Add 50  $\mu$ l of TMB-substrate to each well (also A1) and incubate for 30 min ( $\pm$  2 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.
- 5.10. Stop the reaction by adding 100  $\mu$ l of stop solution to each well (also A1). Positive samples turn yellow.

**Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.**

**The reading should be done within 15 min after adding the Stop Solution!**

#### **5.B. TABLE FOR THE WORKING PROCEDURE**

	Blank (A1)	Negative control	Positive control	Sample
Negative control	-	50 $\mu$ l	-	-
Positive control	-	-	50	-
Sample	-	-	-	50 $\mu$ l
Incubate for 60 min at 37 °C, wash 3 x with 200 $\mu$ l wash buffer				
Antigen-conjugate mixture	-	50/60 $\mu$ l *)	50/60 $\mu$ l *)	50/60 $\mu$ l *)
Incubate for 60 min at 37 °C, wash 3 x with 200 $\mu$ l wash buffer				
TMB-substrate	50 $\mu$ l	50 $\mu$ l	50 $\mu$ l	50
Incubate for 30 minutes at 37 °C in the dark				
Stop solution	100 $\mu$ l	100 $\mu$ l	100 $\mu$ l	100 $\mu$ l
Photometric reading at 450 nm (ref. 620-650 nm)				

\*) manual/automatic procedure (see 5.6.)

#### 6.A. CALCULATION OF RESULTS (VALIDITY)

- \* Read OD values at 450 nm (reference wavelength 620-650 nm).
- \* Subtract the OD value of the blank (well A1) from all other OD values.
- \* The mean OD value of the **negative Control** has to be **< 0.100**.  
The mean OD value of the **positive Control** has to be **> 0.800**.
- \* **Cut-off = mean OD value of the negative control + 0.320**
- \* **Grey zone = Cut-off ± 10 %**

#### 6.B. INTERPRETATION OF RESULTS/LIMITATIONS OF THE METHOD

- \* Samples with OD values below the lower limit of the grey zone are reported as **NEGATIVE**.
- \* Samples with OD values within the grey zone are reported as **EQUIVOCAL**.  
  
These should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- \* Samples with OD values exceeding the upper limit of the grey zone are reported as **POSITIVE**.
- \* The results should always be interpreted in connection with clinical data and additional diagnostic parameters.
- \* In rare cases ANA-positive specimen can cause false-positive results.
- \* Very high concentrations of lipids may lead to higher OD values, which can influence the test result.  
High concentrations of Hemoglobin or Bilirubin do not influence the test results.

**7. PERFORMANCE CHARACTERISTICS**

We determined the following performance characteristics during the diagnostic evaluation.

**7.A. SENSITIVITY AND SPECIFICITY**

A total of 466 serum specimen including healthy blood donors, patients with suspected acute rubella virus infections, with IgM antibodies to other viruses (CMV, HSV, VZV, EBV, Parvo B19), as well as sera containing ANA and rheumatoid factor were tested with the Rubella-IgM-ELA test PCS medac in comparison to a highly specific ELISA that is used routinely in laboratories. The results obtained are shown in the table below.

	Reference test	
	negative	positive
Rubella-IgM-ELA Test PKS medac	360	1
	3	102

Sensitivity = 99.03 %  
 Specificity = 99.14 %

Positive predictive value: 97.14 %  
 Negative predictive value: 99.72 %

**7.B. PRECISION**

Sample	Intra-assay variation				Sample	Inter-assay variation			
	mean OD	SD	CV (%)	n		mean OD	SD	CV (%)	n
NC	0.025	0.004	17	18	NC	0.024	0.005	19	11
PC	1.235	0.031	3	21	PC	0.962	0.057	6	11
N <sup>0</sup> 1	1.400	0.036	3	21	N <sup>0</sup> 4	0.922	0.054	6	11
N <sup>0</sup> 2	0.930	0.026	3	21	N <sup>0</sup> 5	0.230	0.008	3	11
N <sup>0</sup> 3	0.488	0.019	4	21	N <sup>0</sup> 6	1.159	0.066	6	11
					N <sup>0</sup> 7	0.426	0.035	8	11

NC = Negative Control; PC = Positive Control

## **GENERAL HANDLING ADVICES**

- \* Do not exchange the vials and their screw caps in order to avoid cross contamination.
- \* The reagents have to be sealed immediately after use in order to avoid evaporation and microbial contamination.
- \* After use, the reagents have to be stored as indicated to guarantee the shelf life.
- \* After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

## **HEALTH AND SAFETY INFORMATION**

- \* The local occupational safety and health regulations have to be regarded.
- \* Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

## **DISPOSAL CONSIDERATIONS**

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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