cHSP60-IgG-ELISA medac

Deutsch/English/Français

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cHSP60-IgG-ELISA medac

Recombinant enzyme immunoassay for the quantitative detection of
IgG antibodies to chlamydial heat shock protein 60

Cat. no.: 435

FOR IN VITRO DIAGNOSTIC USE ONLY

INTRODUCTION

Heat shock proteins (HSPs) are highly conserved cellular stress proteins; their amino acid composition has not changed very much during evolution and they are expressed in both procaryotes and eucaryotes with pronounced homology. Thus, between men and mice HSP homologies of 99.9 % were found, between men and bacteria of about 60 %, between men and chlamydiae of about 50 %, and between the chlamydia species of > 95 %.

The heat shock protein superfamily includes various families which are classified according to their molecular weight: HSP10, HSP25 (alpha B-crystallin), HSP27, HSP60, HSP70, HSP90, HSP110. The members of the HSP families are divided into two main groups: the constitutive heat shock proteins (HSC) and the inducible heat shock proteins (HSP).

Under physiological conditions the HSCs are permanently produced; they are the intracellular chaperones on molecular level which are responsible for well-regulated courses of anabolism and catabolism.

The HSPs (stress proteins) are rapidly and distinctly synthesized by the cell in response to various physical and chemical stress stimuli in order to take on increased functions with regard to cell protection.

Microbial HSPs (mHSPs) are dominant antigens and in humans extremely immunogenic. During an infection the microbes strongly increase their mHSP synthesis in order to protect themselves from the host's immunological defense mechanisms. In this connection mHSP60 becomes one of the most predominant bacterial proteins. The immunity which follows such primary infections is normally restricted to microbe-specific epitopes of the mHSP60 molecule.

The chlamydial heat shock protein 60 (cHSP60) is one of the most intensively investigated mHSPs. Asymptomatic, not recognized or insufficiently treated chlamydial infections may lead to chronic persistence

of the pathogens. At this stage cHSP60 is continuously overexpressed, i.e., the human immune system permanently has to deal with this foreign protein. Cellular and humoral immune responses are induced. Furthermore, the human HSP60 (hHSP60) is also overexpressed in order to protect the human cells from the chlamydial attacks. The antibodies that are directed towards the conserved epitopes of the cHSP60 may cross-react with the conserved regions of the hHSP60 which finally may lead to autoimmune responses.

An immunity to the conserved epitopes of cHSP60 and, in parallel, a cross-immunity to hHSP60 with resulting possible autoimmunity cannot be recognized without a specific, standardized detection system.


A variety of sequelae resulting from chlamydia infections have been linked with cHSP60-induced autoimmune responses: adverse outcome of pregnancy, preterm labour, habitual and spontaneous abortions, ectopic pregnancy, tubal factor infertility (TFI), and *in vitro* fertilization (IVF) failure.

Other autoimmune diseases, such as reactive arthritis, distinct forms of asthma, as well as arteriosclerosis and coronary heart diseases as possible autoimmune diseases due to sensitization to HSP60, to date, partly have been discussed controversially. The different results, among others, may be put down to the use of not standardized in-house tests, which do not allow direct comparison. A standardized, reproducible, objective test system is mandatory in order to draw clear conclusions from comparable test results in connection with the various clinical pictures.


The cHSP60-IgG-ELISA medac uses recombinant heat shock protein 60 from ***Chlamydia trachomatis*** as antigen. It is suitable for the detection of IgG antibodies to HSP60 from ***Chlamydia trachomatis***. In spite of the high homology between the chlamydia species of > 95 % on protein level this assay seems to detect predominantly antibodies which are directed towards HSP60 from ***Chlamydia trachomatis***.

The cHSP60-IgG-ELISA medac should be used for the diagnosis of the tubal status in the forefront of the *in vitro* fertilization. In combination with the ***Chlamydia trachomatis***-specific MOMP serology from medac it represents an additional marker for the detection of ***Chlamydia trachomatis***-induced tubal damages.


TEST PRINCIPLE


The plate is coated with recombinant heat shock protein 60 from *Chlamydia trachomatis* (cHSP60).

The antibodies from the specimen that are directed towards cHSP60 are bound to the antigen.

Peroxidase-conjugated anti-human IgG antibodies bind to the IgG antibodies (P = peroxidase).

Incubation with TMB-substrate (*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

Advantages of the test

- ☞ Recombinant cHSP60 as antigen.
- ☞ The breakable microwell strips permit efficient use of the test.
- ☞ Suitable for automation on open ELISA devices.

KIT CONTENTS

Cat.-no.: 435

1.

MTP

Microplate: 12 x 8 wells (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with recombinant heat shock protein 60 from *Chlamydia trachomatis* and FCS, ready to use.
2.

CONTROL	-
----------------	----------

Negative control: 1 vial with 1.5 ml, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycine sulfate.
3.

CONTROL	+
----------------	----------

Positive control: 1 vial with 1.5 ml, human serum, ready to use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
4.

WB

Wash buffer: 1 bottle with 100 ml, PBS/Tween (10x), pH 7.2 - 7.4, contains ProClin™ 300.
5.

BAC-DIL

Sample diluent: 1 bottle with 110 ml, PBS/Tween/NBCS, pH 7.0 - 7.2, ready to use, stained blue, contains ProClin™ 300.
6.

CON

Conjugate: 3 vials with 4.5 ml each, goat anti-human IgG, HRP-conjugated, ready to use, stained green, contains BSA, phenol, ProClin™ 300 and gentamycine sulfate.
7.

TMB

TMB-substrate: 1 vial with 10 ml, ready to use.
8.

STOP

Stop solution: 2 vials with 11 ml each, 0.5 M sulfuric acid (H₂SO₄), ready to use.

1. STORAGE AND STABILITY

MATERIAL/REAGENT	STATE	STORAGE	STABILITY
Test kit	unopened	2 - 8 °C	until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	12 weeks
Controls	opened	2 - 8 °C	12 weeks
Wash buffer	diluted	2 - 8 °C	12 weeks
Sample diluent	opened	2 - 8 °C	12 weeks
Conjugate	opened	2 - 8 °C	12 weeks
TMB-substrate	opened	2 - 8 °C	12 weeks
Stop solution	opened	2 - 8 °C	until expiry date

Do not use the reagents after the expiry date.

2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- 2.1. Water for injection (H₂O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

3. PREPARATION OF THE REAGENTS

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

3.1. Microplate

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated under point 1.

3.2. Wash buffer

Mix one volume of wash buffer (10x) with nine volumes of water for injection (e.g., 50 ml wash buffer (10x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

Crystals in the wash buffer (10x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.

Do not mix test specific reagents (microplate, controls, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.

Reagents from other manufacturers must not be used in general.

Valid and reproducible results are only obtained if the test procedure is precisely followed.

4. SPECIMEN

4.1. The test is suitable for serum samples.

4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.

4.3. Sera have to be diluted 1:50 with sample diluent.

5.A. TEST PROCEDURE

5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

Microplate wells are ready to use and do not have to be prewashed.

5.2. Pipette 50 µl of the sample diluent into the well A1 as blank (see 6.A.), and 50 µl of the negative control (in duplicate), positive control, and the diluted patient samples.

If necessary, the microplate wells can be kept in a humid chamber up to 30 min at RT before proceeding.

5.3. Incubate the microplate wells for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.

5.4. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

Do not allow the wells to dry out! Proceed immediately!

5.5. Add conjugate (coloured green) to each well.

50 µl of conjugate have to be pipetted into the wells if the test is done manually.

Please note:

When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.

The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.

5.6. Incubate again for 60 min (\pm 5 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil.

5.7. After incubation wash microplate wells again (see 5.4.).

5.8. Add 50 µl of TMB-substrate to each well and incubate for 30 min (\pm 2 min) at 37 °C (\pm 1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well. Positive samples turn yellow.

Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.

The reading should be done within 15 min after adding the stop solution!

5.B. TABLE FOR THE TEST PROCEDURE

	Blank (A1)	Negative control	Positive control	Sample
Sample diluent	50 µl	-	-	-
Negative control	-	50 µl	-	-
Positive control	-	-	50 µl	-
Sample	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
Conjugate	50/60 µl *)	50/60 µl *)	50/60 µl *)	50/60 µl *)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer.				
TMB-substrate	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark.				
Stop solution	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)				

*) manual/automatic procedure (see 5.5.)

6.A. CALCULATION OF RESULTS (VALIDITY)

- * Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- * Subtract the OD value of the blank (well A1) from all other OD values.
- * The OD value of the blank has to be **< 0.100**.
- * The mean OD value of the **negative control** has to be **< 0.100**.
- * The OD value of the **positive control** has to be **> 0.800**.
- * **Cut off = mean OD value of the negative control + 0.350**
- * **Grey zone = Cut off ± 10 %**

Repeat the run if the results do not meet the specification!

6.B. INTERPRETATION OF THE RESULTS:

6.B.1. QUALITATIVE

Result	Valuation
OD < Grey zone	negative
OD Cut-off +/- 10 %	equivocal
OD > Grey zone	positive

6.B.2. SEMI-QUANTITATIVE

Cut-off Index: $\frac{\text{OD Sample}}{\text{OD Cut-off}}$	Valuation
< 0.9	negative
0.9 - 1.1	equivocal
> 1.1	positive

6.B.3. QUANTITATIVE

The **endtiter** (= last positive or borderline dilution) can be calculated.

First of all the cut-off index has to be calculated.

$$\text{Cut-off index} = \frac{\text{OD sample}}{\text{cut off}}$$

Then the calculation of the endtiter is possible in two ways.

1. $1 : \text{endtiter} = 55.5 \times (\text{OD/cut off})$

The result has to be rounded down to the next lower whole titer step, i.e., to 1:50, 1:100, 1:200, etc.

This calculation is valid for OD values ≤ 2.0 . Samples with OD values > 2.0 should be retested employing a higher starting dilution (e.g., 1:400).

If a higher starting dilution than 1:50 had been used, the above mentioned formula has to be multiplied with the corresponding dilution factor.

Example: Employed sample dilution 1:400, that means an eightfold higher starting dilution.

OD = 1.78
 Cut off = 0.42
1 : endtiter = $55.5 \times (1.78/0.42) \times 8 = 1882 \Rightarrow$ **1:1600**

2. The titer which corresponds to the cut off index can also be determined from the table below.

Cut-off Index	IgG	
	Result	Titer
<0.90	Negative	<1: 50
0.90 - 1.10	Equivocal	
1.11 - 1.80	Positive	1: 50
1.81 - 3.60	Positive	1:100
3.61 - 7.20	Positive	1:200

This calculation is only valid for OD values ≤ 2.0 ; otherwise see above mentioned example.

- * Samples with OD values within the grey zone should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.
- * The results should always be interpreted in connection with anti-**C. trachomatis** IgG and IgA antibodies and with the clinical data of the patients as well as with additional diagnostic parameters.
- * High concentrations of hemoglobin, bilirubin, and of lipids in serum do not have an influence on the results.
- * Cross-reactions with antibodies to other microbial HSP60 cannot be excluded.

6.C. SPECIFIC INTERPRETATION

Possible Results			Interpretation
pELISA		cHSP60	
IgG	IgA	IgG	
+	+	+	High probability of a <i>C. trachomatis</i> -induced TFI
+	-	+	
+	+	-	Indication of a <i>C. trachomatis</i> infection.
+	-	-	Indication of a past <i>C. trachomatis</i> infection.
-	-	+	No serological indication of a <i>C trachomatis</i> infection; cross-reactions with antibodies to HSP60 of other bacteria possible.
-	-	-	No serological indication of a <i>C. trachomatis</i> -infection.

7. PERFORMANCE CHARACTERISTICS

We have determined the following performance characteristics during the validation.

7.A. PREVALENCE

In various cohorts the following cHSP60 IgG antibody prevalence was determined:

Cohort	cHSP60 IgG prevalence
Blood donors	14 % (n=100)
Patients with fertility problems	59 % (n=73)
Patients with rheumatic diseases	64 % (n=44)
Patients with coronary heart diseases	18 % (n=44)

7.B. CLINICAL-DIAGNOSTIC SIGNIFICANCE

Various studies had demonstrated that the cHSP60 serology in combination with the *C. trachomatis*-spezific MOMP serology (pELISA medac) constitutes an additional marker for the detection of *Chlamydia trachomatis*-induced tubal damages.

7.B.1. SENSITIVITIES AND POSITIVE PREDICTIVE VALUES (PPV) WITH REFERENCE TO OCCLUDED TUBES

Authors: Den Hartog J et. al. (2003) Cohorts: IVF patients with (n=59) and without (n=254) distal tubal pathologie (DTP)	Positive antibody results	DTP- present (n=59)	Sens. (%)	PPV (%)
	<i>C. trachomatis</i> - IgG	32/59 (54 %)	54 %	61 %
	cHSP60-IgG	30/59 (51 %)	51 %	44 %
	<i>C. trach.</i> - and cHSP60-IgG	28/59 (47 %)	47 %	68 %
Authors: Clad et al. (2004) Cohorts: IVF patients with patent (n=102) and occluded tubes (n=24)	Positive antibody results	TFI- present (n=24)	Sens. (%)	PPV (%)
	<i>C. trachomatis</i> - IgG	21/24 (87.5 %)	87,5 %	54 %
	cHSP60-IgG	19/24 (79 %)	79 %	39 %
	<i>C. trach.</i> - and cHSP60-IgG	19/24 (79 %)	79 %	63 %
Authors: Surcel et al. (2004) Cohorts: IVF patients with occluded tubes (n=88) and blood donors (n=163)	Positive antibody results	TFI- present (n=88)	Sens. (%)	PPV (%)
	<i>C. trachomatis</i> - IgG	38/88 (43 %)	43 %	63 %
	cHSP60-IgG	52/88 (59 %)	59 %	59 %
	<i>C. trach.</i> - and cHSP60-IgG	30/88 (34 %)	34 %	70 %

7.B.2. SPECIFICITIES AND NEGATIVE PREDICTIVE VALUES (NPV) WITH REFERENCE TO PATENT TUBES

Authors: Den Hartog J et. al. (2003) Cohorts: IVF patients with (n=59) and without (n=254) distal tubal pathologie (DTP)	Negative antibody results	DTP absent (n=254)	Spec. (%)	NPV (%)
	C. trachomatis- IgG	234/254 (92 %)	92 %	90 %
	cHSP60-IgG	216/254 (85 %)	85 %	88 %
	C. trach.- and cHSP60-IgG	241/254 (95 %)	95 %	89 %
Authors: Clad et al. (2004) Cohorts: IVF patients with patent (n=102) and occluded tubes (n=24)	Negative antibody results	TFI- absent (n=102)	Spez. (%)	NPV (%)
	C. trachomatis- IgG	84/102 (82 %)	82 %	97 %
	cHSP60-IgG	72/102 (71 %)	71 %	94 %
	C. trach.- and cHSP60-IgG	91/102 (89 %)	89 %	95 %
Authors: Surcel et al. (2004) Cohorts: IVF patients with occluded tubes (n=88) and blood donors (n=163)	Negative antibody results	TFI- absent (n=163)	Spez. (%)	NPV (%)
	C. trachomatis- IgG	141/163 (87 %)	87 %	74 %
	cHSP60-IgG	127/163 (78 %)	78 %	78 %
	C. trach.- and cHSP60-IgG	150/163 (92 %)	92 %	72 %

7.B. PRECISION

Sample	Intra-assay variation (n=22)			Sample	Interassay variation (n=11)		
	mean OD	SD	CV (%)		mean OD	SD	CV (%)
NC	0.053	0.006	11.0	NC	0.045	0.009	19.1
BC	0.683	0.023	3.3	BC	0.649	0.033	5.0
PC	1.324	0.057	4.3	PC	1.293	0.083	6.4
N° 1	0.067	0.003	4.4	N° 4	0.090	0.021	22.8
N° 2	0.653	0.013	1.9	N° 5	1.047	0.087	8.3
N° 3	1.305	0.047	3.6	N° 6	1.310	0.113	8.7
				N° 7	1.657	0.087	5.2

NC = negative control; BC = weak positive control (not included in the kit);
PC = positive control

GENERAL HANDLING ADVICES

- * To avoid cross contamination do not exchange the vials and their screw caps.
- * The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.
- * After use, the reagents have to be stored as indicated to guarantee the shelf life.
- * After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

HEALTH AND SAFETY INFORMATION

- * The local occupational safety and health regulations have to be regarded.
- * Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

DISPOSAL CONSIDERATIONS

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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LITERATUR/REFERENCES/LITTERATURE

Burian K, Kis Z, Virok D, Endresz V, Prohaszka Z, Duba J, Berencsi K, Boda K, Horvath L, Romics L, Fust G, Gonczol E: Independent and joint effects of antibodies to human heat-shock protein 60 and ***Chlamydia pneumoniae*** infection in the development of coronary atherosclerosis. *Circulation* (2001) 103:1503-1508

Clad A, Petersen EE, Böttcher M: Expanded ***Chlamydia trachomatis*** serology: cHSP60 IgG and its association with tubal occlusion. In: Deák J (ed.) *Proceedings of the European Society for Chlamydia Research*. Pauker Nyomdaipari Kft., Budapest, Hungary (2004) p 86

Den Hartog J, Land J, Stassen F, Slobbe M, Kessels A, Bruggemann C: Predictive value of cHSP60 antibodies in screening for tubal subfertility. 19th Annual Meeting of the European Society for Human Reproduction and Embryo Transfer (EHSRE) June 29-July 2, 2003, Madrid, Spain, xviii183, P-549

Fajac I, Roisman GL, Lacronique J, Polla BS, Dusser DJ: Bronchial gamma delta T-lymphocytes and expression of heat shock proteins in mild asthma. *Eur Respir J* (1997) 10:633-638

Hahn DL, Peeling RW, Dillon E, McDonald R, Saikku P: Serologic markers for ***Chlamydia pneumoniae*** in asthma. *Ann Allergy Asthma Immunol* (2000) 84:227-233

Hoppichler F, Koch T, Dzien A, Gschwandtner G, Lechleitner M: Prognostic value of antibody titre to heat-shock protein 65 on cardiovascular events. *Cardiology* (2000) 94:220-223

Jantos CA, Krombach C, Wuppermann FN, Gardemann A, Bepler S, Asslan H, Hegemann JH, Haberbosch W: Antibody response to the 60-kDa heat-shock protein of ***Chlamydia pneumoniae*** in patients with coronary artery disease. *J Infect Dis* (2000) 181:1700-1705

Kalayoglu MV, Indrawati, Morrison RP, Morrison SG, Yuan Y, Byrne GI: Chlamydial virulence determinants in atherogenesis: the role of chlamydial lipopolysaccharide and heat shock protein 60 in macrophage-lipoprotein interactions. *J Infect Dis* (2000) 181 Suppl 3:S483-489

Kikuta LC, Puolakkainen M, Kuo CC, Campbell LA: Isolation and sequence analysis of the ***Chlamydia pneumoniae*** GroE operon. *Infect Immun* (1991) 59:4665-4669

Larsen B, Birkelund S, Mordhorst CH, Ejstrup L, Andersen LS, Christensen G: The humoral immune response to ***Chlamydia trachomatis*** in patients with acute reactive arthritis. *Br J Rheumatol* (1994) 33:534-540

Lindquist SC, Craig EA: The heat shock proteins. *Annu Rev Genet* (1988) 22:631-636

Mayr C, Richter K, Lilie H, Buchner J: Cpr6 and Cpr7, two closely related HSP90-associated immunophilins from *Saccharomyces cerevisiae*, differ in their functional properties. *J Biol Chem* (2000) 275:34140-34146

Morrison RP, Belland RJ, Lyng K, Caldwell HD: Chlamydial disease pathogenesis. The 57-kD chlamydial hypersensitivity antigen is a stress response protein. *J Exp Med* (1989) 170:1271-1283

Neuer A, Spandorfer SC, Giraldo P, Dieterle S, Rosenwaks Z, Witkin SS: The role of heat shock proteins in reproduction. *Hum Reprod Update* (2000) 6:149-159

Neuer A, Gao Y, Sziller I, Dieterle S: Strategies for extended *Chlamydia trachomatis* serology in infertile patients: A clinical evaluation of *Chlamydia trachomatis* MOMP-, chlamydial 60KD heat shock protein and human HSP60 β serology. In: Deák J (ed.) *Proceedings of the European Society for Chlamydia Research*. Pauker Nyomdaipari Kft., Budapest, Hungary (2004) p. 242

Shinnik TM: Heat shock proteins as antigens of bacterial and parasitic pathogens. *Current Topics Microbiol Immunol* (1991) 167:145-160

Surcel HM, Bloigu A, Öhman H, Halttunen M, Birkelund S, Tiitinen A, Christiansen G, Koskela P, Morrison R, Paavonen J: *C. trachomatis* specific immune responses and infertility. In: Deák J (ed.) *Proceedings of the European Society for Chlamydia Research*. Pauker Nyomdaipari Kft., Budapest, Hungary (2004) p 162

Tishler M, Shoenfeld Y: Anti-heat-shock protein antibodies in rheumatic and autoimmune diseases. *Semin Arthritis Rheum* (1996) 26:558-563

Wang JD, Michelitsch MD, Weissmann JS: GroEL-GroES-mediated protein folding requires an intact central cavity. *Proc Natl Acad Sci USA* (1998) 95:12163-12168

Wick G, Xu Q: Atherosclerosis - an autoimmune disease. *Exp Gerontol* (1999) 34:559-566

Zügel U, Kaufmann SHE: Role of heat shock proteins in protection from and pathogenesis of infectious diseases. *Clin Microbiol Rev* (1999) 12:19-39