**Chlamydia trachomatis-IgA-ELISA *plus* medac**

Deutsch/English/Français

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**Chlamydia trachomatis-IgA-ELISA plus medac**

Enzyme immunoassay for the quantitative detection of IgA antibodies  
to **Chlamydia trachomatis**

Cat. no.: 498-PLUS

FOR IN VITRO DIAGNOSTIC USE ONLY

**INTRODUCTION**

Chlamydiae are gram-negative bacterial pathogens. They have an obligate intracellular life cycle in mucosal surfaces and according to recent findings in endothelial cells and smooth muscle cells. Chlamydiae depend on energy-rich phosphates of their host cells and are therefore called energy parasites.

The genus chlamydia comprises four species: **C. pneumoniae**, **C. psittaci**, **C. pecorum** and **C. trachomatis**. **C. pneumoniae** and **C. trachomatis** are obligate pathogens of humans. **C. psittaci** is pathogenic for humans and a variety of animal species. Up to now **C. pecorum** has been isolated from animals only.

**C. trachomatis** is one of the most frequent agents causing sexually transmitted bacterial infections of the urogenital tract and the eye. During delivery **C. trachomatis** can be transferred to the infant and can cause neonatal conjunctivitis and pneumonia.

**C. trachomatis** is subdivided into 19 serovars causing different clinical manifestations. Serovars A-C are responsible for the endemic trachoma, an ocular disease, which may lead to blindness after years. It affects hundreds of million people in developing countries. Lymphogranuloma venereum, a tropical sexually transmitted disease, is the result of infections with serotypes L<sub>1</sub>-L<sub>3</sub>. In industrial countries serotypes D-K are the cause of the most frequent sexually transmitted oculo-genital infections.

Most **C. trachomatis** infections in women are asymptomatic. This results in the many chronic diseases sustained by the ascended, persisting agents. Subsequent to a peripheral infection in most of the cases cervicitis will develop. Urethritis is less common in women. Between a primary infection and the appearance of chronic sequelae often months or even years will pass.

Asymptomatic *C. trachomatis* infections with mainly chronic course in women may lead to endometritis, adnexitis, periappendicitis, perihepatitis, peritonitis and reactive arthritis. This latter arthritis form is also known in men as so-called post-urethritic arthritis.

Repeated adnexitis may often lead to infertility. In men, ascending infections following asymptomatic urethritis may result in epididymitis and prostatitis. The corresponding inflammatory processes may influence the motility of sperms and the acrosome reaction in a negative way, which results in reduced fertility.

Peripheral *C. trachomatis* pathogens can be detected with high sensitivity and specificity by laboratory methods, such as ELISA and IFT or by nucleic acid amplification techniques like PCR.

In chronic infections serology is method of choice. For the detection of species-specific antibodies the microimmunofluorescence (MIF) has been considered the "gold standard". MIF is laborious, subjective, and requires much experience. This test system is not standardized; the use of different antigens and different cut off criteria for past and recent or current infections creates significant lab-to-lab variations of results.


The ELISA technique is the favoured method in laboratory routine at this time. The quality of the diagnostic information depends as well on the employed antigen as on the method of diagnostic data interpretation. The **Chlamydia trachomatis-IgA-ELISA plus medac** uses a synthetic peptide of a variable domain from an immunodominant region of the MOMP (major outer membrane protein = MOMP).

This highly specific antigen has also been employed in the **Chlamydia trachomatis-IgG-ELISA plus medac**.


The detection of IgA and IgG antibodies to *C. trachomatis* allows a classification of the state of infection.

Moreover, the medac one-point quantification (AU/ml) provides the best preconditions for reproducible results and in this way for measuring of the antibody courses.


## TEST PRINCIPLE


The plate is coated with synthetic *C. trachomatis*-specific peptide.

The *C. trachomatis*-specific antibodies from the specimen are bound to the antigen.

Peroxidase-conjugated anti-human IgA antibodies bind to the IgA antibodies (P = peroxidase).

Incubation with TMB-substrate (\*). The reaction is stopped by the addition of sulfuric acid. The absorption is read photometrically.

### Advantages of the test

- ☞ No false positive results caused by cross-reactivities with other chlamydiae.
- ☞ One-point-quantification, no standard curve needed.
- ☞ The breakable microwell strips permit efficient use of the test.

## KIT CONTENTS

Cat. no.: 498-PLUS

1. 

<b>MTP</b>
------------

Microplate: 12 x 8 wells, purple-coded (with frame and desiccant vacuum sealed in aluminium bag), breakable, U-form, coated with *C. trachomatis*-specific, synthetic peptid and NBCS, ready to use.
2. 

<b>CONTROL</b>	<b>-</b>
----------------	----------

Negative control: 1 vial with 1.5 ml, human serum, ready to use, stained blue, contains NBCS, phenol, ProClin™ 300 and gentamycin sulfate.
3. 

<b>CONTROL</b>	<b>+</b>
----------------	----------

Positive control: 1 vial with 1.5 ml, human serum, ready to use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
4. 

<b>CAL</b>
------------

Calibrator: 1 vial with 1.5 ml, human serum, ready to use, stained blue, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
5. 

<b>WB</b>
-----------

Wash buffer: 1 bottle with 100 ml, PBS/Tween (10x), pH 7.2 - 7.4, contains ProClin™ 300.
6. 

<b>BAC-DIL</b>
----------------

Sample diluent: 1 bottle with 110 ml, PBS/Tween/NBCS, pH 7.0 - 7.2, ready to use, stained blue, contains ProClin™ 300.
7. 

<b>CON</b>
------------

Conjugate: 3 vials with 4.5 ml each, goat anti-human IgA, HRP-conjugated, ready to use, stained yellow, contains BSA, phenol, ProClin™ 300 and gentamycin sulfate.
8. 

<b>TMB</b>
------------

TMB-substrate: 1 vial with 10 ml, ready to use.
9. 

<b>STOP</b>
-------------

Stop solution: 2 vials with 11 ml each, 0.5 M sulfuric acid (H<sub>2</sub>SO<sub>4</sub>), ready to use.

## **1. STORAGE AND STABILITY**

<b>MATERIAL/REAGENT</b>	<b>STATE</b>	<b>STORAGE</b>	<b>STABILITY</b>
Test kit	unopened	2 - 8 °C	until expiry date
Microplate	opened	2 - 8 °C in bag with desiccant	12 weeks
Controls/Calibrator	opened	2 - 8 °C	12 weeks
Wash buffer	diluted	2 - 8 °C	12 weeks
Sample diluent	opened	2 - 8 °C	12 weeks
Conjugate	opened	2 - 8 °C	12 weeks
TMB-substrate	opened	2 - 8 °C	12 weeks
Stop solution	opened	2 - 8 °C	until expiry date

Do not use the reagents after the expiry date.

## **2. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED**

- 2.1. Water for injection (H<sub>2</sub>O redist.). Use of deionised water can disturb the test procedure.
- 2.2. Adjustable micropipettes.
- 2.3. Clean glass or plastic containers for dilution of wash buffer and specimen.
- 2.4. Suitable device for microplate washing (e.g. multistepper or ELISA washer).
- 2.5. Incubator for 37 °C.
- 2.6. Microplate reader with filters for 450 nm and 620 - 650 nm.

## **3. PREPARATION OF THE REAGENTS**

Before starting the test procedure all kit components must be equilibrated to room temperature.

Calculate the number of wells required.

### **3.1. Microplate**

The aluminium bag has to be tightly resealed together with the desiccant after each removal of wells. Storage and stability of the wells are indicated under point 1.

### 3.2. Wash buffer

Mix one volume of wash buffer (10x) with nine volumes of water for injection (e.g. 50 ml wash buffer (10x) with 450 ml water). 10 ml of diluted wash buffer are needed for eight wells.

**Crystals in the wash buffer (10x) have to be dissolved by warming (max. 37 °C) and/or stirring at RT.**

**Do not mix test specific reagents (microplate, controls, calibrator, conjugate) from different kit lots. In contrast to that, sample diluent, wash buffer, TMB-substrate and stop solution are generally exchangeable in all Chlamydia-and Mycoplasma-ELISA medac.**

**Reagents from other manufacturers must not be used in general.**

**Valid and reproducible results are only obtained if the test procedure is precisely followed.**

## **4. SPECIMEN**

4.1. The test is suitable for serum samples.

4.2. Pretreatment of sera, e.g. inactivation, is not necessary. However, they should neither be contaminated with microorganisms nor contain any red blood cells.

4.3. Sera have to be diluted 1:50 with sample diluent.

## **5.A. TEST PROCEDURE**

5.1. Cut the aluminium bag above the zip fastener and take out the required number of microplate wells (see 3.1.).

**Microplate wells are ready to use and do not have to be pre-washed.**

5.2. Pipette 50 µl sample diluent into well A1 as blank (see 6.A.). Add 50 µl each of the negative control, the positive control as well as the diluted samples in single determination, and 50 µl of the calibrator in duplicate to the wells.

**If necessary the microplate wells can be kept up to 30 min at RT before proceeding.**

5.3. Incubate the microplate wells for 60 min (± 5 min) at 37 °C (± 1 °C) in a humid chamber or sealed with incubation cover foil.

5.4. After incubation wash the microplate wells three times with 200 µl wash buffer per well. Pay attention that all wells are filled. After washing tap microplate wells on filter paper.

**Do not allow the wells to dry out! Proceed immediately!**

5.5. Add conjugate (coloured yellow) to each well.

**50 µl of conjugate have to be pipetted into the wells if the test is done manually.**

**Please note:**

**When working with automated devices, 60 µl of conjugate have to be pipetted into each well due to a higher evaporation in the incubation chambers of the devices.**

**The suitability of the test for automated devices was confirmed during the evaluation of the test. Nevertheless we recommend to verify the compatibility of the test with the devices used in the lab.**

5.6. Incubate again for 60 min ( $\pm$  5 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil.

5.7. After incubation wash microplate wells again (see 5.4.).

5.8. Add 50 µl of TMB-substrate to each well and incubate for 30 min ( $\pm$  2 min) at 37 °C ( $\pm$  1 °C) in a humid chamber or sealed with incubation cover foil in the dark. Positive samples turn blue.

5.9. Stop the reaction by adding 100 µl of stop solution to each well. Positive samples turn yellow.

**Clean microplate wells from underneath before the photometric reading and take care that there are no air bubbles in the wells.**

**The reading must be done within 15 min after adding the stop solution!**

### **5.B. TABLE FOR THE TEST PROCEDURE**

	Blank (A1)	Negative control	Positive control	Calibrator	Sample
Sample diluent	50 µl	-	-	-	-
Neg. control	-	50 µl	-	-	-
Pos. control	-	-	50 µl	-	-
Calibrator	-	-	-	50 µl	-
Sample	-	-	-	-	50 µl
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer					
Conjugate	50/60 µl*)	50/60 µl*)	50/60 µl*)	50/60 µl*)	50/60 µl*)
Incubate for 60 min at 37 °C, wash 3 x with 200 µl wash buffer					
TMB-substrate	50 µl	50 µl	50 µl	50 µl	50 µl
Incubate for 30 min at 37 °C in the dark					
Stop solution	100 µl	100 µl	100 µl	100 µl	100 µl
Photometric reading at 450 nm (ref. 620 - 650 nm)					

\*) manual/automatic procedure (see 5.5.)

### **6.A. CALCULATION OF RESULTS (VALIDITY)**

- \* The evaluation is performed using arbitrary units (AU/ml).
- \* Read OD values at 450 nm (reference wavelength 620 - 650 nm).
- \* Subtract the OD value of the blank (well A1) from all other OD values.
- \* Lot-specific data

The lot-specific data sheet provided with the kit contains the following information:

- Lot-specific calibration curve
- Curve parameters a and b
- Nominal OD value of the calibrator
- Lower OD limit of the calibrator
- Nominal concentration range (AU/ml) of the positive control

\* Validity criteria

- The OD value of the **blank** has to be < **0.100**.
- The OD value of the **negative control** has to be < **0.100**.
- The unit value of the **positive control** has to be within the nominal range indicated in the lot-specific data sheet.
- The mean OD value of the **calibrator** has to be above the lower OD limit indicated in the lot-specific data sheet.

**Repeat the run if the results do not meet the specification!**

\* Correction of the results

The measured OD values of the positive control and the samples have to be corrected as follows:

$$OD_{\text{corrected}} = \frac{\text{Nominal OD value of the calibrator}}{\text{Measured OD value of the calibrator}} \times OD_{\text{measured}}$$

\* Quantification of the results

The corresponding concentrations of the corrected OD values in AU/ml can be read from the lot-specific calibration curve (see lot-specific data sheet).

Alternatively, the concentrations can be calculated using the following formula:

$$\text{Concentration [AU/ml]} = b \left( \frac{a}{OD_{\text{corrected}}} - 1 \right)$$

Most of the new ELISA readers allow to program the formula, thus enabling an automated data processing.

The measuring range spans from 22 to 200 AU/ml. Samples below this range have to be interpreted as < 22 AU/ml, those above as > 200 AU/ml. These values must not be extrapolated.

**The cut-off is 25 AU/ml.**

**Grey zone = 22 - 28 AU/ml.**

## 6.B. INTERPRETATION OF RESULTS/LIMITATIONS OF THE METHOD

\* Samples with Unit values below the lower limit of the grey zone are reported as **NEGATIVE**.

\* Samples with Unit values within the grey zone are reported as **EQUIVOCAL**.

These samples should be retested together with a fresh specimen taken 14 days later in order to determine a titer change.

\* Samples with Unit values exceeding the upper limit of the grey zone are reported as **POSITIVE**.

\* The results should always be interpreted in connection with IgG and clinical data and additional diagnostic parameters.

\* High concentrations of hemoglobin, of bilirubin, and of lipids in serum do not have an influence on the results.

**6.C. SPECIFIC IgA-/IgG-INTERPRETATION**

Possible results		Interpretation
IgA AU/ml	IgG AU/ml	
>28	<22	1. Serological indication of early stage of infection or solitary, persisting IgA. <sup>1</sup> . Retest IgA and IgG after 14 days.
>28	>28	2. Serological indication of current infection <sup>2</sup> . Retest IgA and IgG after 14 days.
<22	>28	3. Serological indication of past infection. In case of clinical suspicion retest after 14 days for IgA and IgG antibody courses.
<22	<22	4. No serological indication of current or past infection <sup>3</sup> . In case of clinical suspicion retest after 14 days for IgA and IgG antibodies.

<sup>1</sup> In individual cases solitary IgA antibodies may persist. This immunological phenomenon appears in different bacteriological infections. A clinical relevance is not appraisable.

<sup>2</sup> A current infection may stand for:  
 - chronic infection with persisting pathogens: The unit values of the antibodies remain constant for several weeks.  
 - acute infection: The unit values for the antibodies increase. A twofold increase of the unit value is considered as statistically significant elevated.

<sup>3</sup> In cases of fresh acute *Chlamydia trachomatis* infections the serological antibody results may be negative despite clinical symptoms and positive antigen detection. If a serological confirmation of a positive antigen result or if a follow-up is desired, we recommend to test after 14 days for seroconversion.

**Comment:**

Borderline results may indicate commencing or subsiding infections. Retesting after 14 days is recommended.

## 7. PERFORMANCE CHARACTERISTICS

We determined the following performance characteristics during the diagnostic evaluation

### 7.A. SPECIFICITY AND SENSITIVITY

330 samples were measured in comparison to the Chlamydia trachomatis-IgA-pELISA medac.

The results are shown in the following nine-square table:

#### **C. trachomatis-IgA-pELISA medac**

<b>C. trachomatis IgA-ELISA plus medac</b>		negative	borderline	positive
	negative	243	0	0
	borderline	4	4	2
	positive	3	11	63

Sensitivity = 97%  
Specificity = 97%  
Concordance = 94%

### 7.B. PRECISION

Sample	Intra-assay Variation				Sample	Inter-assay Variation (n = 11)		
	Ø AU/ml	SD	CV (%)	n		Ø AU/ml	SD	CV (%)
PC	92.1	3.3	3.5	22	PC	96.0	4.8	5.0
1	4.2	0.4	9.3	22	5	3.2	0.3	10.9
2	57.8	2.1	3.6	22	6	41.3	2.2	5.3
3	121.1	6.0	5.0	22	7	70.8	3.4	4.8
4	263.6	15.9	6.0	22	8	125.5	9.3	7.4

### **GENERAL HANDLING ADVICES**

- \* To avoid cross-contamination do not exchange the vials and their screw caps.
- \* The reagents have to be sealed immediately after use to avoid evaporation and microbial contamination.

- \* After use, the reagents have to be stored as indicated to guarantee the shelf life.
- \* After use, all components of the testkit should be stored in the original package, in order to avoid mixing up the reagents of other test systems or lots (see also 3.).

#### **HEALTH AND SAFETY INFORMATION**

- \* The local occupational safety and health regulations have to be regarded.
- \* Reagents of human origin have been tested and found to be negative for HBsAg, for antibodies to HIV-1/2 and to HCV. Nevertheless, it is strongly recommended that these materials as well as those of animal origin (see kit contents), should be handled as potentially infectious and used with all necessary precautions.

#### **DISPOSAL CONSIDERATIONS**

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Christiansen, G., Pedersen, L., Clausen, J.D., Birkelund, S.: Cell and molecular biology of Chlamydia. In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, September 11-14, pp. 3-6 (1996).

Fidellers, A.A.A., Land, J.A., Voss, G., Kessels, A.G.H., Severens, J.L.: Cost-effectiveness of Chlamydia antibody tests in subfertile women. Hum. Reprod. 20, 425-432 (2005).

Hoyme, U.B., Spitzbart, H.: Past and current prevalence of ***Chlamydia trachomatis*** in women in Germany. In: Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, September 11-14, p. 391 (1996).

Keck, C., Gerber-Schäfer, C., Clad, A., Wilhelm, C., Breckwoldt, M.: Seminal tract infections: impact on male fertility and treatment options. Hum. Reprod. Update 4, 891-903 (1998).

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Marrazzo, J.M., Johnson, R.E., Green, T.A., Stamm, W.E., Schachter, J., Bolan, G., Hook E.W., Jones, R.B., Martin, D.H., St. Louis, M.E., Black, C.M.: Impact of patient characteristics on performance of nucleic acid amplification tests and DNA probe for detection of ***Chlamydia trachomatis*** in women with genital infections. J. Clin. Microbiol. 43, 577-584 (2005).

Martin, D.H., Nsuami, M., Schachter, J., Hook, E.W., Ferrero, D., Quinn, T.C., Gaydos, C.: Use of multiple nucleic acid amplification tests to define the infected-patient „gold standard“ in clinical trials of new diagnostic tests for ***Chlamydia trachomatis*** infections. J. Clin. Microbiol. 42, 4749-4758 (2004).

Nickel, J.C.: Chronic Prostatitis: An infectious disease? *Infect. Urol.* 13, 31-38 (2000).

Ochsendorf, F.R., Ozdemir, K., Rabenau, H., Fenner, T., Oremek, R., Milbradt, R., Doerr, H.W.: **Chlamydia trachomatis** and male infertility: chlamydia-IgA antibodies in seminal plasma are **C. trachomatis** specific and associated with an inflammatory response. *Eur. Acad. Dermatol. Venereol.* 12, 143-152 (1999).

Paavonen, J.: **Chlamydia trachomatis**: A major cause of mucopurulent cervicitis and pelvic inflammatory disease in women. In: Sexually Transmitted Diseases. Advances in Diagnosis and Treatment. *Curr. Probl. Dermatol. Elsner. P., Eichmann, A. (eds.). Basel, Karger, 24, 110-122 (1996).*

Patton, D., Askienazy-Elbhar, M., Henry-Suchet, J., Campbell, L.A., Cappuccio, A., Tannous, W., Wang, S.P., Kuo, C.C.: Detection of **Chlamydia trachomatis** in fallopian tube tissue in women with post-infectious tubal infertility. *Am. J. Obstet. Gynecol.* 171, 95-101 (1994).

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Petersen, E.E., Clad, A., Mendel, R., Prillwitz, J., Hintz, K.: Prevalence of chlamydial infections in Germany: Screening of asymptomatic women and men by testing first void urine by ligase chain reaction (LCR). In: *Chlamydia Research. Angelika Stary (ed.). Proceedings of the Third Meeting of the European Society for Chlamydia Research, Vienna, Austria, September 11-14, p. 415 (1996).*

Schachter, J.: The intracellular life of Chlamydia. *Curr. Top. Microbiol. Immun.* 138, 109-39 (1988).

Ward, M.E.: The immunobiology and immunopathology of chlamydial infections. *APMIS* 103, 769-796 (1995).

Weström, L.V.: Chlamydia and its effect on reproduction. *J. Brit. Fertil. Soc.* 1, 23-30 (1996).

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